

Mouse Acbp Recombinant Protein

Catalog Number: PROTP31786

Overview

Product Name	Mouse Acbp Recombinant Protein
Description	Mouse Acbp Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-87aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4)

Concentration

1mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Mouse, P31786, 1-87aa; MGSSHHHHHH SSGLVPRGSH MGSMQAQEFK KAAEEVKRLK TQPTDEEMLF IYSHFKQATV GDVNTDRPGL
LDLKGKAKWD SWNKLKGTSK ESAMKTYVEK VDELKKKYGI

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse Acbp Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.