

## Human STI1 Recombinant Protein

Catalog Number: PROTP31948

### Overview

Product Name	Human STI1 Recombinant Protein
Description	Human STI1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-543aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

1mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, P31948, 1-543aa; MEQVNELKEK GNKALSVGNI DDALQCYSEA IKLDPHNHVL YSNRSAAYAK KGDYQKAYED GCKTVDLKPDP  
WGKGYSRKAA ALEFLNRFEE AKRTEEGLK HEANNPQLKE GLQNMEARLA ERKFMNPFNM PNLQKLESD PRTRTLSDP TYRELIEQLR  
NKPSDLGTLK QDPRIMTTLS VLLGVDLGSM DEEEEIATPP PPPPKKETK PEPMEEDLPE NKKQALKEKE LGNDAYKKKD FDTALKHYDK  
AKELDPTNMT YITNQAAYVF EKGDYNKCRE LCEKAIEVGR ENREDYRQIA KAYARIGNSY FKEEKYKDAI HFYNKSLAEH RTPDVLKKCQ  
QAEKILKEQE RLAYINPDLA LEEKNKGNEC FQKGDYPQAM KHYTEAIKRN PKDAKLYSNR AACYTKLLEF QLALKDCEEC IQLEPTFIKG  
YTRKAAALEA MKDYTKAMDV YQKALDLSS CKEAADGYQR CMMAQYNRHD SPEDVKRRAM ADPEVQQIMS DPAMRLILEQ MQKDPQALSE  
HLKNPVIAQK IQKLMDVGLI AIR

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human STI1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.