

Human recombinant CD66b/CEACAM8 protein

Catalog Number: PROTP31997

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human recombinant CD66b/CEACAM8 protein
Description	Human recombinant CD66b/CEACAM8 protein expressed in Baculovirus with His-Tag. Sequence domain: 35-320aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

1mg/ml (determined by Absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, P31997, 35-320aa; QLTIEAVPSN AAEGKEVLLL VHNL PQDPRG YNWKGETVD ANRRIIGYVI SNQQITPGPA YSNRETIYPN ASLLMRNVTR NDTGSYTLQV IKLNL MSEEV TGQFSVHPET PKPSSNNNS NPVEDKDAVA FTCEPETQNT TYLWWVNGQS LPVSPRLQLS NGNRTL TLLS VTRNDVGPYE CEIQNPASAN FSDPVT LNLV YGPDAPTISP SDTYHAGVN LNLSCHAASN PPSQYSWSVN GTFQQYTQKL FIPNITTKNS GSYACHTTNS ATGRNR TTVR MITVSD

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Human recombinant CD66b/CEACAM8 protein](#)

For Research Use Only. Not for use in diagnostic procedures.