

Human HPD Recombinant Protein

Catalog Number: PROTP32754

Overview

Product Name	Human HPD Recombinant Protein
Description	Human HPD Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-393aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 20% glycerol, 50mM NaCl

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P32754, 1-393aa; MTTYSDKGAK PERGRFLHFH SVTFWVGNAK QAASFYCSKM GFEPLAYRGL ETGSREVVSH VIKQ GKIVFV
LSSALNPWNK EMGDHLVKHG DGVKDIAFEV EDCDYIVQKA RERGA KIMRE PWVEQDKFGK VKFAVLQTYG DTHTLVEKM NYIGQFLPGY
EAPAFMDPLL PKLPKCSLEM IDHIVGNQPD QEMVSASEWY LKNLQFHRFW SVDDTQVHTE YSSLRSIVVA NYEESIKMPI NEPAPGKKKS
QIQEYVDYNG GAGVQHIALK TEDIITAIRH LRERGLEFLS VPSTYYKQLR EK LKTAKIKV KENIDALEEL KILVDYDEKG YLLQIFTKPV QDRPTLFLEV
IQRH NHQGFAGNFNSLFKA FEEEQNLRGN LTNMETNGVV PGM

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human HPD Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.