

Human Common beta Chain Recombinant Protein

Catalog Number: PROTP32927

Overview

Product Name	Human Common beta Chain Recombinant Protein
Description	Human Common beta Chain Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 17-443aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.5mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, P32927, 17-443aa; WERSLAGAE TIPLQLRCY NDYTSHTCR WADTQDAQRL VNVTLIRRVN EDLLEPVSCD LSDDMPWSAC
PHPRCVPRRC VIPCQSFVVT DVDYFSFQPD RPLGTRLTVT LTQHVQPPEP RDLQISTDQD HFLTWSVAL GSPQSHWLSP GDLEFEVYK
RLQDSWEDAA ILLSNTSQAT LGPEHLMPS TYVARVRTRL APGSRLSGRP SKWSPEVCWD SQPGDEAQPQ NLECFDGA VLSCSWEVRK
EVASSVSFGL FYKPSPDAGE EECSPVLREG LGS LHTRHHC QIPVPDPATH GQYIVSVQPR RAEKHIKSSV NIQMAPP SLN VTKDGD SYSL
RWETMKMRYE HIDHTEFIQY RKDTATWKDS KTETLQNAHS MALPALEPST RYWARVRVRT SRTGYNGIWS EWSEARSWDT ESVLPMWLEH
HHHHH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Common beta Chain Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.