

Human recombinant ICAM-3/CD50 protein

Catalog Number: PROTP32942-1

Overview

Product Name	Human recombinant ICAM-3/CD50 protein
Description	Human recombinant ICAM-3/CD50 protein expressed in Baculovirus with hlgG-His-Tag. Sequence domain: 30-485aa. Application(s): SDS-PAGE, Bioactivity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hlgG-His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P32942, 30-485aa; QEFLLRV EPQNPVLSAG GSLFVNCSTD CPSSEKIALE TSLSKELVAS GMGWAAFNLS NVTGNSRILC SVYCNGSQIT
GSSNITVYRL PERVELAPLP PWQPVGQNFT LRCQVEDGSP RTSLTVVLLR WEEELSRQPA VEEPAEVTAT VLASRDDHGA PFSCRTELDM
QPQGLGLFVN TSAPRQLRTF VLPVTPPRLV APRFLEVETS WPDVCTLDGL FPASEAQVYL ALGDQMLNAT VMNHGDTLTA TATATARADQ
EGAREIVCNV TLGGERREAR ENLTVFSFLG PIVNLSEPTA HEGSTVTVSC MAGARVQVTL DGVPAAAPGQ PAQLQLNATE SDDGRSFFCS
ATLEVDGEFL HRNSSVQLRV LYGPKIDRAT CPQHLKWKDK TRHVLQCQAR GNPYPELRCL KEGSSREVPV GIPFFVNVTH NGTYQCQASS
SRGKYTLVVV MDIEAGSSH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human recombinant ICAM-3/CD50 protein

For Research Use Only. Not for use in diagnostic procedures.