

## Mouse MMP-2 Recombinant Protein

Catalog Number: PROTP33434

### Overview

Product Name	Mouse MMP-2 Recombinant Protein
Description	Mouse MMP-2 Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 30-662aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

0.25mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 85% by SDS-PAGE

### Amino Acid Sequence

Mouse, P33434, 30-662aa; APSPI IKFPGDVAPK TDKELAVQYL NTFYGCPKES CNLFVLKDTL KKMQKFFGLP QTGDLDQNTI ETMRKPRCGN PDVANYNFFP RPKPKWDKNQI TYRIIGYTPD LDPETVDDAF ARALKVWSDV TPLRFSRIHD GEADIMINFG RWEHGDGYPF DGKDGLLAHA FAPGTGVGGD SHFDDDELWT LGEGQVVRVK YGNADGEYCK FPFLFNGREY SSCTDTGRSD GFLWCSTTYN FEKDGKYGFC PHEALFTMGG NADGQPCKFP FRFQGTSYNS CTTEGRTDGY RWCGTTEDYD RDKKYGFCE TAMSTVGGNS EGAPCVFPFT FLGNKYESCT SAGRNDGKVV CATTNYYDDD RKWGFCDPQG YSLFLVAAHE FGHAMGLEHS QDPGALMAPI YTYTKNFRLS HDDIKGIQEL YGSPDADTD TGTGPTPTLG PVTPEICKQD IVFDGIAQIR GEIFFFKDRF IWRTVTPRDK PTGPLLVATF WPELPEKIDA VYEAPQEEKA VFFAGNEYWV YSASTLERY PKPLTSLGLP PDVQQVDAAF NWSKNKTYI FAGDKFWRYN EVKKKMDPGF PKLIADSWNA IPDNLDVVD LQGGGHSYFF KGAYYKLLEN QSLKSVKFGS IKSDWLGC

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Mouse MMP-2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.