

Human SP-D Recombinant Protein

Catalog Number: PROTP35247

Overview

Product Name	Human SP-D Recombinant Protein
Description	Human SP-D Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 22-375aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.25mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P35247, 22-375aa; EMKTYSH RTMPSACTLV MCSSVESGLP GRDGRDGREG PRGEKGDPLG PGAAGQAGMP GQAGPVGPKG DNGSVGEPGP KGDTGPSGPP GPPGVPGPAG REGPLGKQGN IGPGKPGPK GEAGPKGEVG APGMQGSAGA RGLAGPKGER GVPGERGVPG NTGAAGSAGA MGPQGSPGAR GPPGLKGDKG IPGDKGAKGE SGLPDVASLR QQVEALQGQV QHLQAAFSQY KKVELFPNGQ SVGEKIFKTA GFVKPFTEAQ LLCTQAGGQL ASPRCAAENA ALQQLVVAKN EAAFLSMTDS KTEGKFTYPT GESLVYSNWA PGEPNDDGGS EDCVEIFTNG KWNDRACGEK RLVVCEF

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human SP-D Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.