

Human recombinant TIE1 protein

Catalog Number: PROTP35590-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human recombinant TIE1 protein
Description	Human recombinant TIE1 protein expressed in HEK293 with hlgG-His-Tag. Sequence domain: 22-759aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hlgG-His-Tag
Form	Liquid
Source	HEK293
Formulation	20mM Tris-HCl (pH 8.0) containing 0.1M NaCl, 50% glycerol

Concentration

0.5mg/ml (determined by Absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P35590, 22-759aa; AVDLTLLANL RLTPQRRFFL TCVSGEAGAG RGSDAWGPPL LLEKDDRIVR TPPGPPLRLA RNGSHQVTLR
GFSKPSDLVG VFSCVGGAGA RRTRVIYVHN SPGAHLLPDK VTHTVNKGDV AVLSARVHKE KQTDVIWKS NGSYFYTLTDWH EAQDGRFLLQ
LPNVQPPSSG IYSATYLEAS PLGSAFFRLI VRGCGAGRWG PGCTKECPGC LHGGVCHDHD GECVCPPGFT GTRCEQACRE GRFGQSCQEQ
CPGISGCRGL TFCLPDPYGC SCGSGWGRSQ CQEACAPGHF GADCRLQCQC QNGGTCDRFS GCVCPSGWHG VHCEKSDRIP QILNMASELE

FNLETMPRIN CAAAGNPFPV RGSIELRKP D GTVLLSTKAI VEPEKTTAEF EVPRVLADS GFWECRVSTS GGQDSRRFKV NVKVPPVPLA
APRLLTKQSR QLVVSPLVSF SGDGPSTVR LHYRPQDSTM DWSTIVVDPS ENVTLMLNLRP KTGYSVRVQL SRPGEGGEGA WGPPTLMTTD
CPELLQPWL EGWHVEGTDR LRVSWSLPLV PGPLVGDGFL LRLWDGTRGQ ERRENVSSPQ ARTALLTGLT PGTHYQLDVQ LYHCTLLGPA
SPPAHVLLPP SGPPAPRHLH AQALSDSEIQ LTKWHPEALP GPISKYVVEV QVAGGAGDPL WIDVDRPEET STIIRGLNAS TRYLFMRAS
IQGLGDWSNT VEESTLGNGL QAEGPVQESR AAEEGLDQ

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Human recombinant TIE1 protein](#)

For Research Use Only. Not for use in diagnostic procedures.