

Human BSG Recombinant Protein

Catalog Number: PROTP35613

Overview

| | |
|--------------|---|
| Product Name | Human BSG Recombinant Protein |
| Description | Human BSG Recombinant Protein expressed in Baculovirus with hIgG-His-tag. Sequence domain: 22-207aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method). |
| Size | Starting from 10ug |
| Tag | hIgG-His-Tag |
| Form | Liquid |
| Source | Baculovirus |
| Formulation | Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol |

Concentration

0.5mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P35613, 22-207aa; ADPAAGTVFT TVEDLGSKIL LTCSLNDSAT EVTGHRWLKG GVVLKEDALP GQKTEFKVDS DDQWGEYSCV
FLPEPMGTAN IQLHGPPRVK AVKSSEHINE GETAMLVCKS ESVPVTDWA WYKITDSEDK ALMNGSES RFVSSSQGRSE LH IENLNMEA
DPGQYRCNGT SSKGSDQAI I TLRVRSHLAL EPKSCDKTHT CPPCPAPELL GGPSVFLFPP KPKDTLMISR TPEVTCVVVD VSHEDPEVKF
NWYVDGVEVH NAKTKPREEQ YNSTYRVVSV LTVLHQDWLN GKEYKCKVSN KALPAIEKT ISKAKGQPRE PQVYTLPPSR DELTKNQVSL
TCLVKGFYPS DIAVEWESNG QPENNYKTPP VLVDSDGSFF LYSKLTVDKS RWQQGNVFC SVMHEALHNH YTQKLSLSP GKHHHHHH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human BSG Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.