

Rat Peroxiredoxin 2/Prdx2 Recombinant Protein

Catalog Number: PROTP35704

Overview

Product Name	Rat Peroxiredoxin 2/Prdx2 Recombinant Protein
Description	Rat Peroxiredoxin 2/Prdx2 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-198aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol 1mM DTT

Concentration

1mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Rat, P35704, 1-198aa; MGSSHHHHHH SGLVPRGSH MGSHMASGNA HIGKPAPDFT GTAVVDGAFK EIKLSDYRGK YVVLFFYPLD FTFVCPTEII
AFSDHAEDFR KLGCEVLGVS VDSQFTHLAW INTPRKEGGL GPLNIPLAD VTKLSLQNYG VLKNDEGIAY RGLFIIDAKG VLRQITVNDL
PVGRSVDEAL RLVQAFQYTD EHGEVCPAGW KPGSDTIKPN VDDSKEYFSK HN

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Rat Peroxiredoxin 2/Prdx2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.