

## Rat recombinant TrkA protein

Catalog Number: PROTP35739

### Overview

Product Name	Rat recombinant TrkA protein
Description	Rat recombinant TrkA protein expressed in HEK293 with hlgG-His-Tag. Sequence domain: 35-418aa. Application(s): SDS-PAGE, Bioactivity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hlgG-His-Tag
Form	Liquid
Source	HEK293
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol.

### Concentration

1 mg/ml (determined by Absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Rat, P35739, 35-418aa; SCRETCCPVG PSGLRCTRAG TLNLTGRGRG AGNLTELYVE NQRDLQRLEF EDLQGLGELR SLTIVKSGLR FVAPDAFHFT PRLSHLNLSS NALESLSWKT VQGLSLQDLT LSGNPLHCSC ALLWLQRWEQ EDLCGVYVYQK LQSGSGDQF LPLGHNNSCG VPSVKIQMPN DSVEVGDDVF LQCQVEGQAL QQADWILTEL EGTATMKKSG DLPSLGLTLV NVTSDLNKKV VTCWAENDVG RAEVSVQVSV SFPASVHLGK AVEQHHCIP FSDGQAPAS LRWFFNGSVL NETSFIPTQF LESALTNEM RHGCLRLNQP THVNNNGNYTL LAANPYGQAA ASIMAAFMDN PFEFNPEDPI PVSFSPVDTN STSRDPVEKK DETP

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Rat recombinant TrkA protein

For Research Use Only. Not for use in diagnostic procedures.