

Human Chitinase 3-like Recombinant Protein

Catalog Number: PROTP36222

Overview

Product Name	Human Chitinase 3-like Recombinant Protein
Description	Human Chitinase 3-like Recombinant Protein expressed in CHO Cell with His-Myc-tag. Sequence domain: 1-383aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Myc-Tag
Form	Liquid
Source	CHO Cell
Formulation	Phosphate-Buffered Saline (pH 7.4)

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P36222, 1-383aa; MGVKASQTGF VVLVLLQCCS AYKLVCCYTS WSQYREGDGS CFPDALDRFL CTHIIYSFAN ISNDHIDTWE WNDVTLYGML NTLKNRNP NL KTLSSVGGWN FGSQRFSKIA SNTQSRRTFI KSVPPFLRTH GFDGLDLAWL YPGRRDQKHF TLLIKEMKAE FIKEAQPGKK QLLLSAALSA GKVTIDSSYD IAKISQHLDF ISIMTYDFHG AWRGTTGHHS PLFRGQEDAS PDRFSNTDYA VGYMLRLGAP ASKLVMIPT FGRSFTLASS ETGVGAPISG PGIPGRFTKE AGTLAYEIC DFLRGATVHR ILGQQVPYAT KGNQWVGYDD QESVKSQVQY LKDRQLAGAM VWALDLDDFQ GSFCGQDLRF PLTNAIKDAL AAT

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Chitinase 3-like Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.