

Human PGM1 Recombinant Protein

Catalog Number: PROTP36871

Overview

Product Name	Human PGM1 Recombinant Protein
Description	Human PGM1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-562aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

1mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, P36871, 1-562aa; MGSSHHHHHH SGLVPRGSH MGSMVKIVTV KTQAYQDQKP GTSGLRKRVK VQSSANYAE NFIQSIISTV
EPAQRQEATL VVGDDGRFYM KEAIQLIARI AAANGIGRLV IGQNGILSTP AVSCIIRKIK AIGGIILTAS HNPGGPNGDF GIKFNISNGG PAPEAITDKI
FQISKTIIEY AVCPDLKVDL GVLGKQQFDL ENKFKPFTVE IVDSVEAYAT MLRSIFDFSA LKELLSGPNR LKIRIDAMHG VVGYPYVKKIL CEELGAPANS
AVNCVPLEDF GGHHPDPNLT YAADLVETMK SGEHDFGA AF DGDGDRNMIL GKHGFFVNPS DSVAVIAANI FSIPYFQQTG VRGFARSMPT
SGALDRVASA TKIALYETPT GWKFFGNLMD ASKLSLCGEE SFGTGS DHIR EKDGLWAVLA WLSILATRKQ SVEDILKDHV QKYGRNFFTR
YDYEEVEAEG ANKMMKDLEA LMFDRSFVGK QFSANDKYYT VEKADNFEYS DPVDG SISR N QGLRLIFTDG SRIVFRLSGT GSAGATIRLY
IDSYEKDVAK INQDPQVMLA PLISIALKVS QLQERTGRTA PTVIT

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human PGM1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.