

Human CD96 Recombinant Protein

Catalog Number: PROTP40200

Overview

Product Name	Human CD96 Recombinant Protein
Description	Human CD96 Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 22-503aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.5mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P40200, 22-503aa; VWEKTVNTEE NVYATLGSDV NLTCQTQTVG FFVQMWSKV TNKIDLIQVY HPQYGFYCAV GRPCESLVTF
TETPENGSKW TLHLRNMSCS VSGRYECMLV LYPEGIQTKI YNLLIQTHVT ADEWNSNHTI EIEINQTL EIEINQTL EIEINQTL EIEINQTL
SQNHLSNST LLKDRVKLG T DYRLHLSVQ IFDDGRKFSC HIRVGNPKIL RSSTTVKVFA KPEIPVIVEN NSTDVLEVER FTCLLKNVFP KANITWFIDG
SFLHDEKEGI YITNEERK GK DGFELEKSVL TRVHSNKPAQ SDNLIWCMA LSPVPGNKVW NISSEKITFL LGSEISSTDP PLSVTESTLD TQPSPASSVS
PARYPATSSV TLVDVSALRP NTPQPSNSS MTRGFNYPW TSSGDTKKS VSRISETYS SSPSGAGSTL HDNVFTSTAR AFSEVPTTAN
GSTKTNHVHI TGIVVNKPKD GM

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human CD96 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.