

## Mouse Thrombopoietin Recombinant Protein

Catalog Number: PROTP40226

### Overview

|              |   |
|--------------|---|
| Product Name | Mouse Thrombopoietin Recombinant Protein  |
| Description  | Mouse Thrombopoietin Recombinant Protein expressed in HEK293 with His-tag. Sequence domain: 22-356aa. Application(s): SDS-PAGE, Bioactivity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method). |
| Size         | Starting from 10ug  |
| Tag          | His-Tag   |
| Form         | Liquid  |
| Source       | HEK293  |
| Formulation  | Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol  |

### Concentration

0.5mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Mouse, P40226, 22-356aa; SPVAPACDPR LLNKLRLDSH LLHSRLSQCP DVDPLSIPVL LPAVDFSLGE WKTQTEQSKA QDILGAVSLL  
LEGVMAARGQ LEPSCLSSLL GQLSGQVRLI LGALQGLLGT QLPLQGRRTA HKDPNALFLS LQQLLRGKVR FLLLVEGPTL CVRRTLPTTA  
VPSSTSQLLT LNKFPNRTSG LLETNFSVTA RTAGPGLLSR LQGFRVKITP GQLNQTSRSP VQISGYLNRT HGPVNGTHGL FAGTSLQTLE  
ASDISPGAFN KGSLAFNLQG GLPPSPSLAP DGHTPFPPSP ALPTTHGSP QLHPLFPDPS TTMPNSTAPH PVTMYPHPRN LSQET

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Mouse Thrombopoietin Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.