

## Human STAT3 Recombinant Protein

Catalog Number: PROTP40763

### Overview

Product Name	Human STAT3 Recombinant Protein
Description	Human STAT3 Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 1-770aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 30% glycerol, 1mM DTT

### Concentration

0.25mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, P40763, 1-770aa; MAQWNQLQQL DTRYLEQLHQ LYSDSFPME L RQFLAPWIES QDWAYAASKE SHATLVFHNL LGEIDQQYSR FLQESNVLYQ HNLRRIKQFL QSRYLEK PME IARIVARCLW EESRLLQTAA TAAQGGQAN HPTAAVVTEK QQMLEQLHQD VRKRVQDLEQ KMKVVENLQD DDFDFNYKTLK SQGDMQDLNG NNQSVTRQKM QLEQMLTAL DQMRRSIVSE LAGLLSAMEY VQKTLTDEEL ADWKRRQQIA CIGGPPNICL DRLENWITSL AESQLQTRQQ IKKLEELQK VSYKGDPIVQ HRPMLERIV ELFRNLMKSA FVVERQPCMP MHPDRPLVIK TGVQFTTKVR LLVKFPELNY QLKIKVCIDK DSGDVAALRG SRKFNILGTN TKVMNMEESN NGLSAEFKH LTLREQRCGN GGRANCDASL IVTEELHLIT FETEVYHQGL KIDLETHSLP VVISNICQM PNAWASILWY NMLTNNPKNV NFFT KPIGT WDQVAEVLWSW QFSSTTKRGL SIEQLTTLAE KLLGPGVNYS GCQITWAKFC KENMAGKGF FWWLDNIID LVKKYILALW NEGYIMGFIS KERERAILST KPPGTFLLR F SESSKEGGVT FTWVEKDISG KTQIQSVEPY TKQQLNMSF AEIIMGYKIM DATNILVSPL VYLYPDIPKE EAFGKYCRPE SQEHPEADPG SAAPYLKTKF ICVTPTCSN TIDLPMSPRT LDSLMQFGNN GEGAEPSAGG QFESLTFDME LTSECATSPM

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human STAT3 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.