

## Mouse recombinant FasL (Fas ligand) protein, AF

Catalog Number: PROTP41047

### Overview

Product Name	Mouse recombinant FasL (Fas ligand) protein, AF
Description	Fas Ligand (FasL) is a 17.34 kDa tumor necrosis factor with 152 amino acid residues. FasL is mainly expressed from lymphoid tissue and secreted to blood. Binding to its receptor, TNFRSF6/FAS, leads to induce apoptotic signal into cells. Involved in cytotoxic T-cell-mediated apoptosis, natural killer cell-mediated apoptosis and in T-cell development.
Size	5ug,20ug,100ug
Tag	His Tag (N-term)
Form	Lyophilized
Source	Escherichia coli
Formulation	The protein was lyophilized from a 0.2 um filtered solution containing 1X PBS, pH 8.0. If you have any concerns or special requirements, please confirm with us.

### Storage

Lyophilized protein should be stored at -20°C for 1 year. Upon reconstitution, store at 2°C to 8°C for up to 1 week. Further dilute in a buffer containing a carrier protein or stabilizer (e.g. 0.1% BSA, 10%FBS, 5%HSA or 5% trehalose solution), protein aliquots should be stored at -20°C or -80°C for 3-6 months. Avoid repeated freeze/thaw cycles.

### Purity

>98% as determined by SDS-PAGE.

### Molecular weight

The protein has a calculated MW of 18.27 kDa. The protein migrates as 17-25 kDa under reducing condition (SDS-PAGE analysis).

### Activity

Measure by its ability to induce apoptosis in Jurkat cells. The ED<sub>50</sub> for this effect is <1 ug /mL.

### Endotoxin

<0.1 EU per 1 ug of the protein by the LAL method.

### Amino Acid Sequence

QIANPSTPSEKKEPRSV AHLTG NPHSR SIPLEWEDTYGTALISGVK YKKGGLVINETGLYFVYSKVYFRGQSCNNQPLNHK VYMRNSKYPEDLVL MEEKR  
LNYCTTGQIWAHSSYLGA VFNLT SADHLYVNISQLSLINFEESK TFFGLYKL with polyhistidine tag at the N-terminus.

### Reconstitution

Centrifuge at 3000 rpm for 5 mins before opening. It is recommended to reconstitute the lyophilized protein in sterile H<sub>2</sub>O to a concentration not less than 100 ug/mL and incubate the stock solution at room temperature for at least 20 mins to ensure sufficient re-dissolved. Do Not Vortex! Vigorous shaking may impair the biological activity of the protein.

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse recombinant FasL (Fas ligand) protein, AF

For Research Use Only. Not for use in diagnostic procedures.