

Bcl-6 Protein Human Recombinant Protein

Catalog Number: PROTP41182

Overview

Product Name	Bcl-6 Protein Human Recombinant Protein
Description	Bcl-6 Human Recombinant (aa 1-150) expressed in E.coli, shows a 43kDa band on SDS-PAGE (Including GST tag). The Bcl-6 is purified by proprietary chromatographic techniques.
Size	2ug, 5ug, 10ug
Tag	
Form	Sterile Filtered clear solution.
Source	E. Coli
Formulation	Bcl-6 protein in 50mM Tris-HCl, pH7.5 and 10mM L-glutathione (reduced).

Concentration

Storage

Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.

Purity

!description_after_attributes! ?/description_after_attributes?

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Bcl-6 Protein Human Recombinant Protein