

## Human DAN Recombinant Protein

Catalog Number: PROTP41271

### Overview

Product Name	Human DAN Recombinant Protein
Description	Human DAN Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 18-181aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 10% glycerol, 1mM DTT

### Concentration

1mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, P41271, 18-181aa; MGSSHHHHHH SGLVPRGSH MGSAPPINK LALFPDKSAW CEAKNITQIV GHSGCEAKSI QNRACLGQCF  
SYVNPNTFPQ STESLVHCDS CMPAQSMWEI VTLECPGHEE VPRVDKLVK ILHCSCQACG KEPSHEGLSV YVQGEDGPGS QPGTHPHPHP  
HHPHGGQTPE PEDPPGAPHT EEEGAED

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human DAN Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.