

Human IL-12 R beta 1 Recombinant Protein

Catalog Number: PROTP42701

Overview

Product Name	Human IL-12 R beta 1 Recombinant Protein
Description	Human IL-12 R beta 1 Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 24-545aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 30% glycerol, 1mM EDTA, 0.1mM PMSF.

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, P42701, 24-545aa; CRTSECCFQD PPYPDADSGS ASGPRDLRCY RISSDRYECS WQYEGPTAGV SHFLRCCLSS GRCCYFAAGS
ATRLQFSDQA GVSVLYTVTL WVESWARNQT EKSPEVTLQL YNSVKYEPPL GDKVSKLAG QLRMEWETPD NQVGAEVQFR HRTSPSPWKL
GDCGPQDDDT ESCLCPLMN VAQEFQLRRR QLGSQGSSWS KWSSPVCVPP ENPPQPQVRF SVEQLGQDGR RRLTLKEQPT QLELPEGCQG
LAPGTEVTYR LQLHMLSCPC KAKATRTLHL GKMPYLSGAA YNVAVISSNQ FGPGLNQTWH IPADTHTEPV ALNISVGTNG TTMYPWPARAQ
SMTYCIEWQP VGQDGGGLATC SLTAPQDPDP AGMATYSWSR ESGAMGQEK YITIFASAH PEKLTWSTV LSTYHFGGNA SAAGTPHHVS
VKNHSLDSVS VDWAPLLST CPGVLKEYVV RCRDEDSKQV SEHPVQPTET QVTLGSLRAG VAYTVQVRAD TAWLRGVWSQ PQRFSIEVQV SD

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human IL-12 R beta 1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.