

Human recombinant PRCP protein

Catalog Number: PROTP42785

Overview

| | |
|--------------|---|
| Product Name | Human recombinant PRCP protein |
| Description | Human recombinant PRCP protein expressed in HEK293 with His-Tag. Sequence domain: 22-496aa. Application(s): SDS-PAGE, Enzyme Activity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method). |
| Size | Starting from 10ug |
| Tag | His-Tag |
| Form | Liquid |
| Source | HEK293 |
| Formulation | Phosphate-Buffered Saline (pH 7.4) containing 30% glycerol |

Concentration

0.25mg/ml (determined by Absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, P42785, 22-496aa; LRPALRALGS LHLPTNPTSL PAVAKNYSVL YFQQKVDHFG FNTVKTFNQR YLVADKYWKK NGGSILFYTG
NEGDIWFN NTGFMWDVAE ELKAMLVFAE HRYYGESLPF GDNSFKDSRH LNFLTSEQAL ADFAELIKHL KRTIPGAENQ PVIAIGGSYG
GMLAAWFRMK YPHMVVGALA ASAPIWQFED LVPCGVFMKI VTDFRKSGP HCSESIHRSW DAINRLSNTG SGLQWLTGAL HLCSPLOTSQD
IQHLKDWISE TWVNLAMVDY PYASNFLQPL PAWPIKVVCC YLKNPNVSDS LLLQNIQFAL NVYYNYSQV KCLNISETAT SSLGTLGWSY
QACTEVVMPF CTNGVDDMFE PHSWNLKELS DDCFQQWGVR PRPSWITMY GGKNISHTN IVFSNGELDP WSGGGVTKDI TDTLVAVTIS
EGAHHLDLRT KNALDPMSVL LARSLEVRHM KNWIRDFYDS AGKQ

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human recombinant PRCP protein

For Research Use Only. Not for use in diagnostic procedures.