

Human MCAM Recombinant Protein

Catalog Number: PROTP43121

Overview

Product Name	Human MCAM Recombinant Protein
Description	Human MCAM Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 24-559aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.5mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, P43121, 24-559aa; VPGAEQPAP ELVEVEVGST ALLKCGLSQS QGNLSHVDWF SVHKEKRTL IFRVRQGQGS EPGEYEQRLS LQDRGATLAL TQVTPQDERI FLCQGKRPRS QEYRIQLRVY KAPEEPNIQV NPLGIPVNSK EPEEVATCVG RNGYPIQVI WYKNGRPLKE EKNRVHIQSS QTVESSGLYT LQSILKAQLV KEDKDAQFYC ELNYRLPSGN HMKESREVTV PVFYPTEKVV LEPEVGMMLK EGDRVEIRCL ADGNPPPHFS ISKQNPSTRE AEEETTNDNG VLVLEPARKE HSGRYECQGL DLDTMISLLS EPQELLVNYV SDVRVSPAAP ERQEGSSLT TCEAESSQDL EFQWLREETG QVLERGPVLQ LHLDKREAGG GYRCVASVPS IPGLNRTQLV NVAIFGPPWM AFKERKVVVK ENMVLNLSCE ASGHPRPTIS WNVNGTASEQ DQDPQRVLST LNLVLTPELL ETGVECTASN DLGKNTSILF LELVNLTLT PDSNTTGLS TSTASPHTRA NSTSTERKLP EPESRGAAAL EHHHHHHH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human MCAM Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.