

Human Biotinidase/BTD Recombinant Protein

Catalog Number: PROTP43251

Overview

Product Name	Human Biotinidase/BTD Recombinant Protein
Description	Human Biotinidase/BTD Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 44-545aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 4ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.25mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, P43251, 44-545aa; AHTGEESVAD HHEAEYYVAA VYEHPSILSL NPLALISRQE ALELMNQNLDIYEQQVMTAA QKDVQIIVFP
EDGIHGFnFT RTSIYPFLDF MPSPQVVRWN PCLEPHRFND TEVLQRLSCM AIRGDMFLVA NLGKPECHS SDPRCPKDGR YQFNTNVVFS
NNGTLVDRYR KHNLYFEAAF DVPLKVDLIT FDTPFAGRFG IFTCFDILFF DPAIRVLRDY KVKHVVPYPTA WMNQLPLLA IEIQKAFVA
FGINVLAANV HHPVLGMTGS GIHTPLESFW YHDMENPKSH LIIAQVAKNP VGLIGAENAT GETDPSHSKF LKILSGDPYC EKDAQEVHCD
EATKWNVNAP PTFHSEMMYD NFTLVVWGK EGYLHVCSNG LCCYLLYERP TLSKELYALG VFDGLHTVHG TYYIQVCALV RCGGLGFDTG
GQEITEATGI FEFHLWGNFS TSYIFPLFLT SGMTEVPDQ LGWENDHYFL RKSRLSSGLV TAALYGRLYE RDLEHHHHHH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Biotinidase/BTD Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.