

Human recombinant OX40/TNFRSF4 protein

Catalog Number: PROTP43489-2

Overview

Product Name	Human recombinant OX40/TNFRSF4 protein
Description	Human recombinant OX40/TNFRSF4 protein expressed in Baculovirus with hIgG-His-Tag. Sequence domain: 29-214aa. Application(s): SDS-PAGE, Bioactivity. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hIgG-His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.5mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P43489, 29-214aa; LHCVGDYTPS NDRCCHECRP GNGMVSRCR SQNTVCRPCG PGFYNDVVSS KPCKPCTWCN LRSGERKQL
CTATQDTVCR CRAGTQPLDS YKPGVDCAPC PPGHFSPGDN QACKPWTNCT LAGKHTLQPA SNSSDAICED RDPPATQPQE TQGPPARPIT
VQPTEAWPRT SQGPSTRPVE VPGGRA

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human recombinant OX40/TNFRSF4 protein

For Research Use Only. Not for use in diagnostic procedures.