

Human PBEF/Visfatin Recombinant Protein

Catalog Number: PROTP43490

Overview

Product Name	Human PBEF/Visfatin Recombinant Protein
Description	Human PBEF/Visfatin Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-491aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing, 10% glycerol, 0.1mM DTT

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, P43490, 1-491aa; MNPAEAEFN ILLATDSYKV THYKQYPPNT SKVYSYFECR EKKTENSCLR KVKYEETVFY GLQYILNKYL KGKVVTKEDI
QEAKDVYKEH FQDDVFNEKG WNYILEKYDG HLPKIKAVP EGFVIPRGNV LFTVENTDPE CYWLTNWIET ILVQSWYPIT VATNSREQKK
ILAKYLLETS GNLDGLEYKL HDFGYRGVSS QETAGIGASA HLVNFKGTDV VAGLALIKKY YGKDPVPGY SVPAAEHSTI TAWGKDHEKD
AFEHIVTQFS SVPVSVSDS YDIYNACEKI WGEDLRHLIV SRSTQAPLII RPDSGNPLDT VLKVLKILGK KFPVTENSKG YKLLPPYLRV IQGDGVDINT
LQEIIEGMKQ KMWSIENIAF GSGGGLLQKL TRDLLNCSFK CSYVVTNGLG INVFKDPVAD PNKRKSKGRL SLHRTPAGNF VTLEEGKGD
EEYGQDLLHT VFKNGKVTKS YSFDEIRKNA QLNIELEAAH H

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human PBEF/Visfatin Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.