

## Tbp7 (PSMC4) (NM\_006503) Human Recombinant Protein

Catalog Number: PROTP43686

### Overview

Product Name	Tbp7 (PSMC4) (NM_006503) Human Recombinant Protein
Description	Recombinant protein of human proteasome (prosome, macropain) 26S subunit, ATPase, 4 (PSMC4), transcript variant 1
Source	HEK293T
Formulation	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

### Storage

Store at -80°C.

### Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

## Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Tbp7 (PSMC4) (NM\_006503) Human Recombinant Protein