

Human MMP-13 Recombinant Protein

Catalog Number: PROTP45452

Overview

Product Name	Human MMP-13 Recombinant Protein
Description	Human MMP-13 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 104-471aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 0.15M NaCl, 1mM DTT

Concentration

0.5mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence


Human, P45452, 104-471aa; YNVFPRT LKWSKMNLTY RIVNYTPDMT HSEVEKAFKK AFKVWSDVTP LNFTRLHDGI ADIMISFGIK EHGDFYPPFDG
PSGLLAHAFP PGPNYGGDAH FDDDETWTSS SKGYNLFLVA AHEFGHSLGL DHSKDPGALM FPIYTYTGKS HFMLPDDDVQ GIQSLYGP GD
EDPNPKHPKT PDKCDPSLSL DAITSRGET MIFKDRFFWR LHPQQVDAEL FLTKSFWPEL PNRIDAAAYEH PSHDLIFIFR GRKFWALNGY
DILEGYPKKI SELGLPKEVK KISAAVFHED TGKTLFSGN QVWRYDDTNH IMDKDYPRLI EEDFPGIGDK VDAVYEKNGY IYFFNGPIQF
EYSIWSNRIV RVMPANSILW C

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human MMP-13 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.