

## USP5 (NM\_003481) Human Recombinant Protein

Catalog Number: PROTP45974

### Overview

Product Name	USP5 (NM_003481) Human Recombinant Protein
Description	Recombinant protein of human ubiquitin specific peptidase 5 (isopeptidase T) (USP5), transcript variant 2
Source	HEK293T
Formulation	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol

### Storage

Store at -80°C.

### Purity

> 80% as determined by SDS-PAGE and Coomassie blue staining

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



USP5 (NM\_003481) Human Recombinant Protein