

Human CrkL Recombinant Protein

Catalog Number: PROTP46109

Overview

| | |
|--------------|---|
| Product Name | Human CrkL Recombinant Protein |
| Description | Human CrkL Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-303aa. Application(s): SDS-PAGE. |
| Size | Starting from 20ug |
| Tag | His-Tag |
| Form | Liquid |
| Source | E. coli |
| Formulation | 20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 2mM DTT, 0.1M NaCl. |

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P46109, 1-303aa; MSSARFDSSD RSAWYMGVPS RQEAQTRLQG QRHGMFLVRD SSTCPGDYVL SVSENSRVSH YIINSLPNRR
FKIGDQEFDH LPALLEFYKI HYLDTTTLIE PAPRYSPPM GSVSAPNLPT AEDNLEYVRT LYDFPGND AE DLPFKKGEIL VII EKPEEQW
WSARNKGRV GMIPVPYVEK LVRSSPHGKH GNRNSNSYGI PEPAHAYAQP QTTTLP PAVS GSPGAAITPL PSTQNGPVFA KAIQKRVPCA
YDKTALALEV GDIVKVTRMN INGQWEGEVN GRKGLFPFTH VKIFDPQNP D ENE

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human CrkL Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.