

Mouse Glutathione S-transferase P 2 Recombinant Protein

Catalog Number: PROTP46425

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Mouse Glutathione S-transferase P 2 Recombinant Protein
Description	Mouse Glutathione S-transferase P 2 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-210aa. Application(s): SDS-PAGE.
Size	Starting from 4ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

1mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Mouse, P46425, 1-210aa; MGSSHHHHHH SSGLVPRGSH MGSMPPTYIV YFSPGRCEA MRMLLADQGQ SWKEEVVTID TWMQGLLKPT CLYGQLPKFE DGDLTLYQSN AILRHLGRSL GLYGKNQREA AQVDMVNDGV EDLRGKYGTM IYRNYENGKN DYVKALPGHL KPFETLLSQN QGGKAFIVGD QISFADYNLL DLLLIHQVLA PGCLDNFPLL SAYVARLSAR PKIKAFLSSP EHVNRPINGN GKQ

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Mouse Glutathione S-transferase P 2 Recombinant Protein](#)

For Research Use Only. Not for use in diagnostic procedures.