

## Human Glycogenin 1 Recombinant Protein

Catalog Number: PROTP46976

### Overview

Product Name	Human Glycogenin 1 Recombinant Protein
Description	Human Glycogenin 1 Recombinant Protein expressed in E. coli with His-T7-tag. Sequence domain: 1-333aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-T7-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 10% glycerol

### Concentration

1mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, P46976, 1-333aa; MHHHHHHMAS MTGGQQMGRD LYDDDDKDRW GSMTDQAFVT LTTNDAYAKG ALVLGSSLKQ HRTTRRLVVL  
ATPQVSDSMR KVLETVFDEV IMVDVLDSGD SAHLTLMKRP ELGVTTLKLH CWSLTQYSKC VFMDADTLVL ANIDDLFDRE ELSAAPDPGW  
PDCFNSGVFV YQPSVETYNQ LLHLASEQGS FDGGDQGILN TFFSSWATTD IRKHLPIYN LSSISIYSYL PAFKVFVGASA KVVHFLGRVK  
PWNYTYDPKT KSVKSEAHDP NMTHPEFLIL WWNIFFTNLV PLLQQFGLVK DTCSYVNVED VSGAISHLSL GEIPAMAQPF VSSEERKERW  
EQGQADYMGA DSFDNIKRKL DTYLQ

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Human Glycogenin 1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.