

## Mouse Aldh2 Recombinant Protein

Catalog Number: PROTP47738

### Overview

Product Name	Mouse Aldh2 Recombinant Protein
Description	Mouse Aldh2 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 20-519aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 20% glycerol, 1mM DTT

### Concentration

0.5mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 85% by SDS-PAGE

### Amino Acid Sequence

Mouse, P47738, 20-519aa; MGSSHHHHHH SGLVPRGSH MGSSAAATSA VPAPNHQPEV FCNQIFINNE WHDAVSRKTF PTVNPSTGEV ICQVAEGNKE DVDKAVKAAR AAFQLGSPWR RMDASDRGRL LYRLADLIER DRTYLAALET LDNGKPYVIS YLVDLDMVLK CLRYAGWAD KYHGKTIPID GDFFSYRHE PVGVCGQIIP WNFPLLMQAW KLGALATGN VVVMKVAEQT PLTALYVANL IKEAGFPPGV VNIVPGFGPT AGAAIASHEG VDKVAFTGST EVGHLIQVAA GSSNLKRVT L ELGGKSPNII MSDADMDWAV EQAHFALFFN QGQCCCAGSR TFVQENVYDE FVERSVARAK SRVVG NPFDS RTEQGPQVDE TQFKKILGYI KSGQQEGAKL LCGGAAADR GYFIQPTVFG DVKDGMTIAK EEIFGPMQI LKFKTIEEVV GRANDSKYGL AAVFTKDL D KANYLSQALQ AGTVWINCYD VFQAQSPFGG YKMSGSGREL GEYGLQAYTE VKT VTKVPQ KNS

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Mouse Aldh2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.