

Mouse recombinant OX40/TNFRSF4 protein

Catalog Number: PROTP47741

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Mouse recombinant OX40/TNFRSF4 protein
Description	Mouse recombinant OX40/TNFRSF4 protein expressed in HEK293 with hIgG-His-Tag. Sequence domain: 20-211aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	hIgG-His-Tag
Form	Liquid
Source	HEK293
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.25mg/ml (determined by Absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Mouse, P47741, 20-211aa; VTARRL NCVKHTYPSG HKCCRECQPG HGMVSRCDHT RDTLCHPCET GFYNEAVNYD TCKQCTQCNH
RSGSELKQNC TPTQDTVCRG RPTGTPRQDS GYKLGVDVCP CPPGHFSPGN NQACKPWTNC TSLGKQTRHP ASDSLDAVCE DRSLLATLLW
ETQRPTRPT TVQSTTVWPR TSELPSPPTL VTPEGP

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Mouse recombinant OX40/TNFRSF4 protein](#)

For Research Use Only. Not for use in diagnostic procedures.