

## Human ATP5O Recombinant Protein

Catalog Number: PROTP48047

### Overview

|              |   |
|--------------|---|
| Product Name | Human ATP5O Recombinant Protein   |
| Description  | Human ATP5O Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 24-213aa. Application(s): SDS-PAGE. |
| Size         | Starting from 20ug  |
| Tag          | His-Tag   |
| Form         | Liquid  |
| Source       | E. coli   |
| Formulation  | 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 40% glycerol, 0.2M NaCl   |

### Concentration

1mg/ml (determined by Bradford assay)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.

### Purity

> 95% by SDS-PAGE

### Amino Acid Sequence

Human, P48047, 24-213aa; FAKLVRPPV QVYGIEGRYA TALYSAASKQ NKLEQVEKEL LRVAQILKEP KVAASVLNPY VKRSIKVKSL NDITAKERFS  
PLTTNLIINLL AENGRLSNTQ GVVSAFSTMM SVHRGEVPCT VTSASPLEEA TLSELKTVLK SFLSQGQVLK LEAKTDPSIL GGMIVRIGEK  
YVDMSVKTKI QKLGRAMREI V [?/description\\_after\\_attributes?](#)

## Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human ATP5O Recombinant Protein