

Mouse Hpd Recombinant Protein

Catalog Number: PROTP49429

Overview

| | |
|--------------|--|
| Product Name | Mouse Hpd Recombinant Protein |
| Description | Mouse Hpd Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-393aa. Application(s): SDS-PAGE. |
| Size | Starting from 10ug |
| Tag | His-Tag |
| Form | Liquid |
| Source | E. coli |
| Formulation | 20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol |

Concentration

0.5mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Mouse, P49429, 1-393aa; MTTYNNK GPKPERGRFL HFHSVTFWVG NAKQAASFYC NKMGEPLAY RGLTGSREV VSHVIKQGKI VFVLCALNP WNKEMGDHLV KHGDGVKDIA FEVEDCDHIV QKARERGAKE VREPWVEQDK FGKVKFAVLQ TYGDTTHTLV EKINYTGRL PGFEPTYKD TLLPKLPRCN LEIIDHIVGN QPDQEMQSAS EWYLNQLFH RFWSVDDTQV HTEYSSLSI VVTNYEESIK MPINEPAPGR KKSQIQEYVD YNGGAGVQHI ALKTEDIITA IRHLRERGTE FLAAPSYYK LLRENLSAK IQVKESMDVL EELHILVDYD EKGYLLQIFT KPMQDRPTLF LEVIQRHNNHQ GFGAGNFNSL FKAEEEEQAL RGNLTDLEPN GVRSGM

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse Hpd Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.