

Human SARS Recombinant Protein

Catalog Number: PROTP49591-1

Overview

Product Name	Human SARS Recombinant Protein
Description	Human SARS Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-514aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 10% glycerol, 100mM NaCl

Concentration

0.5mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Human, P49591, 1-514aa; MVLDDL FRVDKGGDPA LIRETQEKRF KDPGLVDQLV KADSEWRRCR FRADNLNKLK NLCSKTIGEK MKKKEPVGDD
ESVPENVLSF DDLTADALAN LKVSQIKKVR LLIDEAILKC DAERIKLEAE RFENLREIGN LLHPSVPISN DEDVDNKVER IWGDCTVRKK
YSHVDLVMV DGFEGEKGA VAGSRGYFLK GVLVFLEQAL IQYALRTLGS RGYIPIYTPF FMRKEVMQEV AQLSQFDEEL YKVIKKGSEK
SDDNSYDEKY LIATSEQPIA ALHRDEWLRP EDLPIKYAGL STCFRQEVGS HGRDTRGIFR VHQFEKIEQF VYSSPHDNKS WEMFEEMITT
AEEFYQSLGI PYHIVNIVSG SLNHAASKKL DLEAWFPGSG AFRELVSCSN CTDYQARRLR IRYGQTKKMM DKVEFVHMLN ATMCATTRTI
CAILENYQTE KGITVPEKLG EFMPPGLQEL IPFVKPAPIE QEPSKKQKKQ HEGSKKAAA RDVTLENRLQ NMEVTD A

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human SARS Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.