

## Human SARS Recombinant Protein

Catalog Number: PROTP49591-1

### Overview

Product Name	Human SARS Recombinant Protein
Description	Human SARS Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-514aa. Application(s): SDS-PAGE.
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 10% glycerol, 100mM NaCl

### Storage

Can be stored at +2°C to +8°C for 1 week. For long term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles.

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

MVLDDL FRVDKGGDPA LIRETQEKRF KDPGLVDQLV KADSEWRRRCR FRADNLNKLK NLCSKTIGEK MKKKEPVGDD ESPPENVLSF  
DDLADALAN LKVSQIKKVR LLIDEAILKC DAERIKLEAE RFENLREIGN LLHPSVPISN DEDVDNKVER IWGDCTVRKK YSHVDLVVMV  
DGFEGEKGA VAGSRGYFLK GVLVFLEQAL IQYALRTLGS RGYIPIYTPF FMRKEVMQEV AQLSQFDEEL YKVIKKGSEK SDDNSYDEKY  
LIATSEQPIA ALHRDEWLRP EDLPIKYAGL STCFRQEVGS HGRDTRGIFR VHQFEKIEQF VYSSPHDNKS WEMFEEMITT AEEFYQSLGI  
PYHIVNIVSG SLNHAASKKL DLEAWFPGSG AFRELVSCSN CTDYQARRLR IRYGQTKKMM DKVEFVHMLN ATMCATTRTI CAILENYQTE  
KGITVPEKLEK EFMPPGLQEL IPFVKPAPIE QEPSKKQKKQ HEGSKKKAAA RDVTLENRLQ NMEVTD

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human SARS Recombinant Protein