

Human ARL4D Recombinant Protein

Catalog Number: PROTP49703

Overview

Product Name	Human ARL4D Recombinant Protein
Description	Human ARL4D Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-201aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 2mM DTT, 20% glycerol, 200mM NaCl

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P49703, 1-201aa; MGSSHHHHHH SGLVPRGSH MGNHITEMAP TASSFLPHFQ ALHVVVIGLD SAGKTSLLYR LKFKEFVQSV
PTKGFNTEKI RVPLGSGRGI TFQVWDVGGQ EKLRLWRSY TRRTDGLVFV VDAAEAERLE EAKVELHRIS RASDNQGVVP LVLANKQDQP
GALSAAEVEK RLAVRELAAL TLTHVQGCSA VDGLGLQOGL ERLYEMILKR KKAARGGKKR R

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human ARL4D Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.