

Human recombinant Thrombospondin 5/COMP protein

Catalog Number: PROTP49747-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human recombinant Thrombospondin 5/COMP protein
Description	Human recombinant Thrombospondin 5/COMP protein expressed in HEK293 with His-Tag. Sequence domain: 21-757aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	HEK293
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.5mg/ml (determined by Absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, P49747, 21-757aa; QGQSPL GSDLGPQMLR ELQETNAALQ DVRELLRQQV REITFLKNTV MECDACGMQQ SVRTGLPSVR PLLHCAPGFC FPGVACIQTE SGARCGPCPA GFTGNNGSHCT DVNECNAHPC FPRVRCINTS PGFRCEACPP GYSGPTHQGV GLAFKANKQ VCTDINECET GQHNCVPNSV CINTRGSFQC GPCQPGFVGD QASGCQRRRAQ RFCPDGSPSE CHEHADCVLE RDGSRSCVCA VGWAGNGILC GRDLDGFP DEKLRCPERQ CRKDNCVTVP NSGQEDVDRD GIGDACPDPA DGDGVPNEKD NCPLVRNPDQ RNTDEDKWGD ACDNCRSQKN DDQKDTDQDG

RGDACDDDDID GDRIRNQADN CPRVPNSDQK DSDGDGIGDA CDNCPQKSNP DQADVDHDFV GDACDSQDQ DGDGHQDSRD NCPTVPNSAQ
EDSDHDGQGD ACDDDDDDNDG VPDSRDNCRL VPNPGQEDAD RDGVGDVCQD DFDADKVVDK IDVCPENAEV TLTDFAFQT VVLDPEGDAQ
IDPNWVVLNQ GREIVQTMNS DPGLAVGYTA FNGVDFEGTF HVNTVTDDDY AGFIFGYQDS SSFYVVMWKQ MEQTYWQANP FRAVAEPIQ
LKAVKSSTGP GEQLRNALWH TGDTSQVRL LWKDPRNVGW KDKKSYRWFL QHRPQVGYIR VRFYEGPELV ADSNVLDTT MRGGRLGVFC
FSQENIIWAN LRYRCNDTIP EDYETHQLRQ A

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Human recombinant Thrombospondin 5/COMP protein](#)

For Research Use Only. Not for use in diagnostic procedures.