

## Human GSK-3 beta Recombinant Protein

Catalog Number: PROTP49841

### Overview

Product Name	Human GSK-3 beta Recombinant Protein
Description	Human GSK-3 beta Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 1-420aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 0.5mM PMSF and 30% glycerol

### Concentration

0.25mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 90% by SDS-PAGE

### Amino Acid Sequence

Human, P49841, 1-420aa; MSGRPRTTSF AESCKPVQQP SAFGSMKVSR DKDGSKVTTV VATPGQGPD R PQEVSYTDTK VIGNSGFGVV YQAKLCDSGE LVAIKKVLQD KRFKNRELQI MRKLDHCNIV RLRYFFYSSG EKKDEVYLN L VLDYVPETVY RVARHYSRAK QTLPIYVKL YMYQLFRSLA YIHSFGICHR DIKPQNLLLD PDTAVLKLCD FGS AKQLVRG EPNVSYICSR YYRAPELIFG ATDYTSSIDV WSAGCVLAEL LLGQPIFPGD SGVDQLVEII KVLGTPTR EQ IREMNP NYTE FKFPQIK AHP WTKVFRP RTP PEAIALCSRL LEYTP TARLT PLEACAHSFF DELRDPNVKL PNGRDTPALF NFFTQELSSN PPLATILIPP HARIQAAAST PTNATAASDA NTGDRGQTNN AASASASNST

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human GSK-3 beta Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.