

## Mouse Arsa Recombinant Protein

Catalog Number: PROTP50428

### Overview

Product Name	Mouse Arsa Recombinant Protein
Description	Mouse Arsa Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 18-506aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

### Concentration

0.5mg/ml (determined by absorbance at 280nm)

### Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

### Purity

> 95% by SDS-PAGE

### Amino Acid Sequence

Mouse, P50428, 18-506aa; ADPSPPNILL IFADDLGYGD LGSYGHPSST TPNLDQLAEG GLRFTDFYVP VSLCTPSRAA LLTGRLPVRS GMYPGVLGPS SQGGLPLEEV TLAEVLAARG YLTGMAGKWH LGVGPEGAFL PPHQGFHRFL GIPYSHDQGP CQNLTCFPPD IPCKGGCDQG LVPIPLLNL TVEAQPWLP GLEARYVSFS RDLMADAQRQ GRPFFLYYAS HHTHY PQFSG QSFTKRSGRG PFGDSLME LD GAVGALMTTV GDLGLEETL VIFTADNGPE LMRMSNGGCS GLLRCGKGT FEGGVREPAL VYWPGHITPG VTHELASSLD LLPTLAALTG APLPNVTLDG VDISPLLGT GKSPRKS VFF YPPYPDEI HG VFAVRNGKYK AHFFTQGSAH SDTTS DPACH AANRLTAHEP PLYDLSQDP GENYNVLESI EGV SPEALQA LKHIQLLKAQ YDAAMTFGPS QIAKGEDPAL QICCP SCTP HPV CCHCPGS QSHHHHHH

## Usage

---

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

## Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Mouse Arsa Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.