

Mouse Lumican Recombinant Protein

Catalog Number: PROTP51885

Overview

Product Name	Mouse Lumican Recombinant Protein
Description	Mouse Lumican Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 19-338aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 20% glycerol

Concentration

0.5mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

Mouse, P51885, 19-338aa; QYYDYDIPLF MYGQISPNC APECNCPHSYP TAMYCDDLKL KSVPMVPPGI KYLYLRNNQI DHIDEKAFEN
VTDLQWLILD LLENSKIK GKVFSK LKQL KKLHINYNNL TESVG LPLPKS LQDLQLTNNK ISKLG SFDGL VNLTFIYLQH NQLKEDAVSA SLKGLKSLEY
LDLSFNQMSK LPAGLPTSLL TLYLDNNKIS NIPDEYFKRF TGLQYLRLSH NELADSGVPG NSFNISLLE LDLSYNKLKS IPTVNENLEN YYLEVNELEK
FDVKSFCKIL GPLSYSKIKH LRLDGNPLTQ SSLPPDMYEC LRVANEITVN

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Mouse Lumican Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.