

Human Importin alpha 2/KPNA2 Recombinant Protein

Catalog Number: PROTP52292

Overview

Product Name	Human Importin alpha 2/KPNA2 Recombinant Protein
Description	Human Importin alpha 2/KPNA2 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-529aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 20% glycerol, 1mM DTT

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P52292, 1-529aa; MSTNEN ANTPAARLHR FKNKGKDSTE MRRRRIEVNV ELRKAKKDDQ MLKRRNVSSF PDDATSPLQE NRNNQGTVNW SVDDIVKGIN SSNVENQLQA TQAARKLLSR EKQPPIDNII RAGLIPKFVS FLGRTDCSPI QFESAWALTN IASGTSEQTK VVVDGGAIPA FISLLASPHA HISEQAVWAL GNIAGDGSVF RDLVIKYGAV DLLLALLAVP DMSSLACGYL RNLTWLSNL CRNKNPAPPI DAVEQILPTL VRLHHDDPE VLADTCWAIS YLTDGPNERI GMVVKTGVVP QLVKLLGASE LPIVTPALRA IGNIVTGTDE QTQVVIDAGA LAVFPSLLTN PKTNIQKEAT WTMSNITAGR QDQIQQVNVH GLVPFLVSVL SKADFKTQKE AVWAVTNYTS GGTVEQIVYL VHCGIIEPLM NLLTAKDTKI ILVILDAISN IFQAAENLGE TEKLSIMIEE CGGLDKIEAL QNHENESVYK ASLSLIEKYF SVEEEDQNV VPETTSEGYT FQVQDGAPGT FNF

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Importin alpha 2/KPNA2 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.