

e.coli IdhA Recombinant Protein

Catalog Number: PROTP52643

Overview

Product Name	e.coli IdhA Recombinant Protein
Description	E. coli IdhA Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-329aa. Application(s): SDS-PAGE.
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 10% glycerol, 100mM NaCl

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 95% by SDS-PAGE

Amino Acid Sequence

E.coli, P52643, 1-329aa; MGSSHHHHHH SSGLVPRGSH MGSHEMGLAVY STKQYDKKYL QQVNESFGFE LEFFDFLLTE KTAKTANGCE
AVCIFVNDG SRPVLEELKK HGVKYIALRC AGFNNVDLDA AKELGLKVVR VPAYDPEAVA EHAIGMMMTL NRRIHRAVQR TRDANFSLEG
LTGFTMYGKT AGVIGTGKIG VAMLRILKGF GMRLAFDPY PSAAALELGV EYVDLPTLFS ESDVISLHCP LTPENYHLLN EAAFEQMKNQ
VMIVNTRGA LIDSQAAIEA LKNQKIGSLG MDVYENERDL FFEDKSNDVI QDDVFRRLSA CHNVLFTHGQ AFLTAEALTS ISQTTLQNLN
NLEKGETCPN ELV

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



e.coli IdhA Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.