

Human ARFIP1 Recombinant Protein

Catalog Number: PROTP53367

Overview

Product Name	Human ARFIP1 Recombinant Protein
Description	Human ARFIP1 Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-341aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	Phosphate-Buffered Saline (pH 7.4) 30% glycerol, 1mM DTT

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P53367, 1-341aa; MGSSHHHHHH SGLVPRGSH MGSMASQSPK NSAAEIPVTS NGEVDDSDREH SFNRDLKHSL PSGLGLSETQ
ITSHGFDNTK EGVEAGAFQ GGQRTQTKSG PVILADEIKN PAMEKLELVR KWSLNTYKCT RQIISEKLGR GSRTVDLELE AQIDILRDNK
KKYENILKLA QTLSTQLFQM VHTQRQLGDA FADLSLSLE LHEEFGYNAD TQKLLAKNGE TLLGAINFFI ASVNTLVNKT IEDTLMTVKQ
YESARIEYDA YRTDLEELNL GPRDANTLPK IEQSQHLFQA HKEKYDKMRN DVSVKLKFLE ENKVKVLHNQ LVLFHNAIAA YFAGNQKQLE
QTLKQFHIKL KTPGVDAPSW LEEQ

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human ARFIP1 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.