

Human recombinant JNK3/MAPK10 protein

Catalog Number: PROTP53779-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Human recombinant JNK3/MAPK10 protein
Description	Human recombinant JNK3/MAPK10 protein expressed in Baculovirus with His-Tag. Sequence domain: 1-464aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 30% glycerol, 1mM DTT, 0.1mM PMSF

Concentration

0.25mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 85% by SDS-PAGE

Amino Acid Sequence

Human, P53779, 1-464aa; MSLHFLYYCS EPTLDVKIAF CQGFDKQVDV SYIAKHYNMS KSKVDNQFYS VEVGDSTFTV LKRYQNLKPI
GSGAQGIVCA AYDAVLDRNV AIKKLSRPFQ NQTHAKRAYR ELVLMKCVNH KNIISLLNVF TPQKTLEEFQ DVYLMELMD ANLCQVIQME
LDHERMSYLL YQMLCGIKHL HSAGIIHRDL KPSNIVVKSD CTLKILD FGL ARTAGTSFMM TPYVVTRYR APEVILGMGY KENVDIWSVG
CIMGEMVRHK ILFPGRDYID QWNKVIEQLG TPCPEFMKKL QPTVRNYVEN RPKYAGLTFP KLFPDSL FPA DSEHNK LKAS QARDLLSKML

VIDPAKRISV DDALQHPYIN VWYDPAEVEA PPPQIYDKQL DEREHTIEEW KELIYKEVMN SEEKTKNGVV KGQPSPSGAA VNSSESLLPS
SSVNDISSMS TDQTLASDTD SSLEASAGPL GCCR

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

[Human recombinant JNK3/MAPK10 protein](#)

For Research Use Only. Not for use in diagnostic procedures.