

Human Ataxin 3 Recombinant Protein

Catalog Number: PROTP54252

Overview

Product Name	Human Ataxin 3 Recombinant Protein
Description	Human Ataxin 3 Recombinant Protein expressed in E. coli (non-tagged). Sequence domain: 1-370aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 10ug
Tag	Non-Tagged
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 7.5) containing 2mM DTT, 50mM NaCl, 10% glycerol

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P54252, 1-370aa; MESIFHEKQE GSLCAQHCLN NLLQGEYFSP VELSSIAHQL DEEERMMAE GGVTSEDYRT FLQQPSGNMD
DSGFFSIQVI SNALKVWGLE LILFNSPEYQ RLRIDPINER SFICNYKEHW FTVRKLKQW FNLNSLLTGP ELISDTYLAL FLAQLQQEGY SIFVVKGDLP
DCEADQLQM IRVQQMHRPK LIGEELAQLK EQRVHKTDL RVLEANDGSG MLDEDEEDLQ RALALSRQEI DMEDEEADLR RAIQLSMQGS
SRNISQDMTQ TSGTNLTSEE LKRREAYFE KQQKQKQQQQ QQQQQQQQQQ QQQQGDLSGQ SShPCERPAT SSGALGSDLG DAMSEEDMLQ
AAVTMSLETV RNDLKTEGKK

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human Ataxin 3 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.