

Human YARS Recombinant Protein

Catalog Number: PROTP54577

Overview

Product Name	Human YARS Recombinant Protein
Description	Human YARS Recombinant Protein expressed in E. coli with His-tag. Sequence domain: 1-528aa. Application(s): SDS-PAGE.
Size	Starting from 10ug
Tag	His-Tag
Form	Liquid
Source	E. coli
Formulation	20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 10% glycerol, 0.1M NaCl

Concentration

1mg/ml (determined by Bradford assay)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P54577, 1-528aa; MGDAPSPEEK LHLITRNLQE VLGEEKLKEI LKERELKIYW GTATTGKPHV AYFVPMSKIA DFLKAGCEVT ILFADLHAYL DNMKAPWELL ELRVSYENV IKAMLESIGV PLEKLFKFIG TDYQLSKEYT LDVYRLSSVV TQHDSKKAGA EVVKQVEHPL LSGLLYPGLQ ALDEEYLKVD AQFGGIDQRK IFTFAEKYLP ALGYSKRVHL MNPMVPGLTG SKMSSEEEES KIDLLDRKED VKKKLKKAFK EPGNVENNGV LSFIKHVLFK LKSEFVILRD EKWGGNKTYT AYVDLEKDFK AEVVHPGDLK NSVEVALNKL LDPIREKFNT PALKKLASAA YPDPKQKPM AKGPAKNSEP EEVIPSRLDI RVGKIITVEK HPDADSLYVE KIDVGEAEPR TVVSGLVQFV PKEELQDRLV VVLCNLKPKQ MRGVESQGM L LCASIEGINR QVEPLDPPAG SAPGEHVFK GYEKGQDDE LKPKKKVFEK LQADFKISE CIAQWKQTNF MTKLGSISCK SLKGGNIS

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human YARS Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.