

Human EPHB4 Recombinant Protein

Catalog Number: PROTP54760

Overview

Product Name	Human EPHB4 Recombinant Protein
Description	Human EPHB4 Recombinant Protein expressed in Baculovirus with His-tag. Sequence domain: 16-539aa. Application(s): SDS-PAGE. Endotoxin: < 1 EU per 1ug of protein (determined by LAL method).
Size	Starting from 20ug
Tag	His-Tag
Form	Liquid
Source	Baculovirus
Formulation	Phosphate-Buffered Saline (pH 7.4) containing 10% glycerol

Concentration

0.25mg/ml (determined by absorbance at 280nm)

Storage

Can be stored at 2°C to 8°C for 1 week. For long-term storage, aliquot and store at -20°C to -80°C. Avoid repeated freeze-thaw cycles. (Ships with gel ice.)

Purity

> 90% by SDS-PAGE

Amino Acid Sequence

Human, P54760, 16-539aa; LEETLLNTKL ETADLKWVTF PQVDGQWEEL SGLDEEQHSV RTYEVCDVQR APGQAHWLRT GWVPRRGAVH VYATLRF TML ECLSLPRAGRSCKETFTVFY YESDADTATA LTPAWMENPY IKVDTVAAEH LTRKRPGAEA TGKVNKTLR LGPLSKAGFY LAFQDQGACM ALLSLHLFYK KCAQLTVNLT RFPETVPREL VVPVAGSCVW DAVPAPGPSP SLYCREDGQW AEQPVTGCSC APGFEEAEGN TKCRACAQGT FKPLSGEGSC QPCPANSN TIGSAVCQR VGYFRARTDP RGAPCTPPS APRSVVSR LN GSSLHLEWSA PLESGGREDL TYALRCRECR PGGSCAPCGG DLTFDPGPRD LVEPWVVVRG LRPDFTYTFE VTALNGVSSL ATGPVPFEPV NVTTDREVPP AVSDIRVTRS SPSSLSLAWA VPRAPSGAVL DYEVKYHEKG AEGPSSVRF L KTSENRAELR GLKRGASYLV QVRARSEAGY GPFGEHHSQ TQLDESEGWR EQLAVEHHHH HH

Usage

Boster's products are furnished for LABORATORY RESEARCH USE ONLY. The product may not be used as drugs, agricultural or pesticidal products, food additives or household chemicals.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Human EPHB4 Recombinant Protein

For Research Use Only. Not for use in diagnostic procedures.